

Fully Customizable Combinatorial Libraries for Focused Screening of Variants

Tavneet Gill, Daniel Ramirez, Caitlin Hoeber, Su Maw, Lucy Jia Xu, Tony Cox, Elian Lee


Abstract

Combinatorial DNA libraries are a collection of unique synthetic DNA sequences that are mutated simultaneously across targeted regions of a protein and assembled in precise ratios and combinations. These libraries can be cloned and propagated into microorganisms to perform high throughput screening for the identification of synergistic and beneficial mutations, that ultimately enhance organismal function. Combinatorial libraries are widely used in a variety of applications such as protein optimization (binding, stability, activity, and expression) via directed evolution or in vitro antibody affinity maturation via phage display libraries.



Scientists are often screening redundant sequences in their libraries, as a result of unavoidable side products generated by alternative fabrication technologies on the market. Utilizing Twist Bioscience's silicon DNA synthesis platform, precise mutations are individually synthesized and introduced via high quality oligo pools with low 1:2000 error rate and 90% within 5x of the Mean uniformity with no print limit. Coupled with our molecular biology expertise in gene synthesis and DNA assembly, Twist Bioscience can build full-length, highly diverse and precise libraries to maximize your screening capabilities.

Combinatorial DNA Libraries offered at Twist Bioscience are robust and uniquely customizable. By offering multiple combinatorial DNA library varieties, in-silico designs, codon optimization, NGS sequencing and numerous other customizations, we can help accelerate your design-test-build cycle.

Workflow




Domain Customization




Combinatorial Library Domains

(A) Wild Type Biased Ratio library. Traditional Method for reducing number of overall combinatorial substitutions from wild type. An alternative would be the precision library where users can select number of single/double/triple substitutions away from the wild type. **(B) Uniform Variant Library.** Equal Representation of each desired AA Change in each position with uniform distribution at all positions **(C) High Diversity Uniform Ratio Library.** Comprehensive sampling of combinatorial diversity in a domain **(D) Length Variation.** Exploration of variations in length at defined ratios and various amino acid profiles

In-Silico Removal of Liabilities




Leveraging Precision Synthesis



New Products

Spread-Out Low Diversity (SOLD) Library
A tool designed to enable the exploration of the sequence space encoded by mutating multiple positions simultaneously scattered across the parental sequence to generate a diverse library of proteins for enzyme evolution, enzyme engineering, enzyme stabilization, protein engineering and protein stability.



Cloned Oligo Pools
Cloned oligo pools enable the high-throughput screen of variants introduced from an oligo pool. Applications include sgRNA libraries, promoter and barcode libraries or single domain CDR3 libraries. Cloned oligo pools are amplicon sequenced to ensure tight uniformity and coverage with minimal dropouts. **(A/B) Post Cloning Distribution and GC Representation.** **(C) Read Counts before and after Cloning**

