

Development of an optimized workflow for sensitive variant detection in FFPE-damaged samples

Yvain Desplat, Sean Tighe, Katsuya Colón, Michael Bocek, Mathieu Chauleau, Tong Liu, Lydia Bonar, Yufeng Qian, Tina Han, Owen Smith, Elian Lee, Esteban Toro, Siyuan Chen
Twist Bioscience, South San Francisco, CA, USA



Abstract

Background/Objectives: Next-Generation Sequencing (NGS) of formalin-fixed and paraffin-embedded (FFPE) samples allows for detailed characterization of cancer from routinely collected clinical specimens. However, the preservation process often leads to significant damage and degradation of genetic material, complicating library preparation. Furthermore, because sample availability is typically limited, there is a high demand for robust workflows that can efficiently convert more molecules into sequenceable libraries. Given the high variability in damage and yield across FFPE samples, generating high-quality NGS data for variant and indel identification requires a unified, high-efficiency workflow.

Methods: Presented here is a workflow leveraging the Twist TrueAmp library preparation kit and target enrichment workflow to maximize the conversion of damaged and low mass FFPE samples for NGS.

Results: We demonstrate library conversion of severely damaged (DIN<2) FFPE with as little as 5ng input mass. Compared to other workflows, an increase in library complexity is observed. Importantly, the enzymatic fragmentation of input DNA can be tuned to provide larger insert sizes that are appropriate for sequencing with longer read length. Coupled with target enrichment, more target sites are captured and the uniformity is significantly improved. These improvements are driven by the incorporation of Twist TrueAmp polymerase, an engineered polymerase with a low error rate capable of amplifying through GC extremes to recover challenging genomic regions.

Conclusion: In summary, our optimized library preparation kit built on engineered enzymes emerges as a valuable asset for the deployment of NGS-based FFPE assays.

TrueAmp Library Preparation Kit Workflow




Figure 1. Method Outline for Twist TrueAmp Library Preparation Kit Library preparation relies on optimized enzymes: a tunable enzymatic fragmentation module, an efficient ligation module utilizing Twist engineered T4 DNA ligase, and a faithful amplification module with Twist TrueAmp Polymerase.




Figure 2. Tunable Enzymatic Fragmentation Fragmentation time titration demonstrates how different insert sizes can be obtained depending on downstream application needs.

