

Industry-Leading Anti-Idiotypic Antibody Discovery Services

Twist Biopharma Solutions delivers high-quality, critical antibody reagents, derived from Abveris' hyperimmune DiversimAb™ mouse platform. Reagent generation and development is a critical component of preclinical studies and timely clinical advancement. With 250+ successful anti-idiotypic discovery projects, we are an experienced biopharmaceuticals partner to de-risk your reagent development efforts. Twist Biopharma Solutions has crafted a tailored anti-idiotypic antibody discovery platform that offers hybridoma and single B-cell workflows to fit diverse assay applications (PK, ADA, and beyond) with rapid, predictable discovery timelines.

Working with Twist Biopharma Solutions

KEY ADVANTAGES


- 250+ projects completed
- Ultra-rapid, predictable timelines
- Flexible deliverables
- Royalty-free and fee-for-service models
- Validated hybridoma and B-cell workflows
- 96% success rate*

ADVANCED SCREENING CAPABILITIES

- Matrix interference assessment
- Ligand blocking/non-blocking screening
- Antibody pair identification
- Informed affinity assessments
- On-cell binding

STATE-OF-THE-ART TECHNOLOGIES

- Bruker Beacon Select™
- Carterra® LSA
- Sartorius
- Fortebio
- Hamilton
- Tecan



ANTI-ID: Our anti-id guarantee promises delivery of anti-idiotypic antibodies capable of specific binding to the target biologic in the presence of human serum.

For Research Use Only. Not for use in diagnostic procedures.


*based on standard screening protocol

Rapid hybridoma

BENEFITS

- Industry-leading screening package
- Highly diverse and highly qualified antibody deliverables
- 10-12 weeks to sequences, 14-16 weeks to purified antibody supply
- Direct access to experienced project management

WORKFLOW




Ultra-rapid B-cell screening

BENEFITS

- Fastest speed for aggressively tight timelines
- Reliable and sensitive screening for precision in candidate selection
- 6-8 weeks to sequences, 10-12 weeks to purified antibody supply
- Direct access to experienced project management

WORKFLOW



LEARN MORE

[www.twistbioscience.com/
twist-biopharma-solutions](http://www.twistbioscience.com/twist-biopharma-solutions)

biopharmasales@twistbioscience.com

TESTIMONIAL

"At Promega, we struggled to acquire a high-affinity, high-specificity antibody against one target of interest. After failing with two previous vendors, we enlisted the help of Twist Biopharma Solutions and their DiversimAb approach, and we soon had a multitude of high-quality clones to choose from. We eventually chose a clone that far exceeded our expectations for performance. We were continually impressed by their high-quality analytical data, excellent communication, and overall collaborative spirit."

—Promega Corporation